




# PATIENT MONITORING SUITE

The **ITTS Patient Monitoring Suite of Simulators** can be used in all phases of patient care for any human (role-playing) or **ANY** patient simulator.

## ITTS Capnograph Simulator



*Designed to be used as a mainstream carbon dioxide monitor that can be connected in-line with a mask or tube.*

*Simulator displays include carbon dioxide partial pressure, respiratory rate, and CO2 waveform.*

## ITTS Pulse Oximeter Simulator



*The ITTS Pulse Oximeter simulator utilizes the widely used finger-style probe for monitoring blood oxygen saturation.*

*Simulator display includes Heart Rate, Respiratory Rate, SpO2, and SpO2 waveforms.*

## ITTS NIBP Monitor Simulator



*Designed to replicate a portable NIBP monitor, this simulator allows medical responders to gather vital signs without instructor intervention. Onboard controls allow various modes of operation. Display includes systolic & diastolic pressures (mmHg) and heart rate.*

EMAIL: [INFO@TOMMANIKIN.COM](mailto:INFO@TOMMANIKIN.COM)

TOLL FREE 1.866.839.3455



Selected as a  
**2015 Hot Product**  
by the Editors of JEMS

# TOMMANIKIN.COM